

Proprietary and Confidential.

DOCUMENT NAME: Meter Interface Control Document

D/N: Rev:

FILE NAME: BG Meter Update: 2018/02/23

Revision History

Date. Done

by. The update. New rev.

1 2011/05/12 Carol Primary Edition. Ver1.0

2 2011/07/12 Carol List TD-4222d not supported comm and Ver1.1

3 2011/12/15 Carol Modification Ver1.2

4 2012/4/26 Carol Modification. Ver1.3

5 2013/3/12 Carol Command 0x26 modification. Ver1.4

6 2013/6/3 Carol Command 0x54 modification. Ver1.5

7 2013/09/27 Carol Command 0x26 modification. Ver1.6

8 2017/03/28 Carol Note for different Bluetooth versions Ver1.7

9 2018/02/23 Carol Modification. Ver1.8

Document Name: BG Meter multiple parameters Update: 2018/02/23

Document Version:1. 8

Proprietary and confidential. 2 / 16

Table Of Contents

1. GENERAL	
1.1. OVERVIEW	
1.2. DOCUMENT OVERVIEW	
1.3. DEFINITIONS , ACRONYMS , AND ABBREVIATIONS	
2. INTERFACE SPECIFICATION	
2.1. FRAME STRUCTURE	
2.2. COMMAND LIST	
2.2.1. [0 X23] READ DEVICE CLOCK TIME (4 BYTES)	
2.2.2. [0 X24] READ DEVICE MODEL	
2.2.3. [0 X25] READ THE STORAGE DATA WITH INDEX , PART 1(TIME)	
2.2.4. [0 X26] READ THE STORAGE DATA WITH INDEX , PART 2(RESULT)	
2.2.5. [0 X27] READ DEVICE SERIAL NUMBER , PART 1.....	
2.2.6. [0 X28] READ DEVICE SERIAL NUMBER , PART 2.....	
2.2.7. [0 X2B] READ STORAGE NUMBER OF DATA	
2.2.8. [0 X33] WRITE SYSTEM CLOCK TIME (4 BYTES)	
2.2.9. [0 X50] TURN OFF THE DEVICE	
2.2.10. [0 X52] CLEAR /D ELETE ALL MEMORY	
2.2.11. [0 X54] NOTIFICATION FOR ENTERING COMMUNICATION MODE	

Document Name: BG Meter multiple parameters Update: 2018/02/23

Document Version:1. 8

Proprietary and confidential. 3 / 16

1. General

1.1. Overview

The document describes the commands which will be used to communicate with blood glucose meter.

Before retrieving data in the device, please make sure you have already established the proper connection channel between the device and the host side. For example, if you use USB cable, be sure that the USB driver has been properly installed; or if you use Bluetooth communication, be sure that the connection has been successfully established.

Generally speaking, the device will be the slave side waiting for the request from the host which is usually the computer, smart phone, or the hub/gateway.

Note for different Bluetooth versions:

For Bluetooth v2 meter: Please connect to meter by SPP (Serial Port Profile). Meter always stands by for both pairing and connection, the PIN for pairing is 111111.

For Bluetooth v3 meter: Please connect to meter by SPP (Serial Port Profile). Meter must be paired before any connection, please find the instruction in meter manual for the pairing process. The meter supports SSP (Simple Secure Pairing), as a result PIN is not required. If host does not support SSP, please use 0000 as PIN if required.

For Bluetooth v4 meter: Please connect to meter by following GATT service:

UUID Base: 1212-efde-1523-785feabcd123

Service UUID: 0x1523

Characteristic: 0x1524 (write/notify)

1.2. Document Overview

For the first time user on the device integration, please read the 2.1 frame structure. If you are familiar with serial port communication, you can review the command list table in section 2.2 directly. The following sub sections are all about the details for each command.

Document Name: BG Meter multiple parameters Update: 2018/02/23

Document Version:1. 8

Proprietary and confidential. 4 / 16

1.3. Definitions, Acronyms, and Abbreviations

GW Gateway / Hub / Computer side

MD Medical Device, ie. blood glucose meter, blood pressure meter or scales, etc.

CMD Command. Usually its a one byte alphanumeric symbol.

ACK Acknowledgement. Usually its a one byte alphanumeric symbol.

NA Not available

Document Name: BG Meter multiple parameters Update: 2018/02/23

Document Version:1. 8

Proprietary and confidential. 5 / 16

2. Interface specification

Block Diagram

The communication interface is based on RS232 serial port communication :

2.1. Frame Structure

Byte 1 2 3 ~ n-2 n-1 n

Name Start CMD Index, Address, or Data Stop Check
Sum

Format HEX

GW MD 51 CMD A3 [1n-1]

GW MD 51 ACK A5 [1n-1]

ie. Check sum means the summation from byte 1 to byte n-1 with 8-bit length

Example : 8-byte frame

Byte 1 2 3 4 5 6 7 8

Name Start CMD Data_0 Data_1 Data_2 Data_3 Stop ChkSum

Format HEX

GW MD 51 CMD Data_0 Data_1 Data_2 Data_3 A3 [1..7]

GW MD 51 ACK Data_0 Data_1 Data_2 Data_3 A5 [1..7]

ie. Check sum means the summation from byte 1 to byte 7 with 8-bit length Baud rate 19200

Data bits 8

Parity No

Start bit 1

Stop bit 1 Medical Device /

Sensor

(MD) Gateway / Hub /

Computer

(GW) TX

RX

Document Name: BG Meter multiple parameters Update: 2018/02/23

Document Version:1. 8

Proprietary and confidential. 6 / 16

2.2. Command List

Message name CMD/ACK Direction Length

(Bytes) Note

Read device clock time 0x23 GW MD 8

0x23 GW MD 8

Read device model 0x24 GW MD 8

0x24 GW MD 8

Read storage data, part 1(time) 0x25 GW MD 8

0x25 GW MD 8

Read storage data , part 2(result) 0x26 GW MD 8

0x26 GW MD 8

Read device serial number, part 1 0x27 GW MD 8

0x27 GW MD 8

Read device serial number, part 2 0x28 GW MD 8

0x28 GW MD 8

Read storage number of data 0x2B GW MD 8

0x2B GW MD 8

Write device clock time 0x33 GW MD 8

0x33 GW MD 8

Turn off the device 0x50 GW MD 8

0x50 GW MD 8

Clear/Delete all memory 0x52 GW MD 8

0x52 GW MD 8

Notification for entering communication

mo de (only support cable transmission) 0x54 GW MD 8

Note: TD-4222d not support command 0x24 and 0x54

Document Name: BG Meter multiple parameters Update: 2018/02/23
Document Version:1. 8

Proprietary and confidential. 7 / 16

2.2.1. [0x23] Read device clock time (4 bytes)

Message Name: Read device clock time (4 bytes)

Message ID: 0x23

Message Description: Read device date and time (to minute)

Src / Des: GW / MD

Length: Request: 8 bytes / Response: 8 bytes

Format

Message name Direction CMD Data_0 Data_1 Data_2 Data_3

Read System Clock time . GW MD 0x23 0x0 0x0 0x0 0x0

Response from device. GW MD 0x23 Day + Month + Year Minute Hour

Field Description

Field Size

(Bytes) Value Description.

CMD 1 0x23 GW commands MD to response device clock time

DATA 4 0x0 NA

Field Size

(Bytes) Value Description.

ACK 1 0x23 MD response device clock time

DATA 2 Data_1 + Data_0 : Day+Month+Year (word) (see Table A)

DATA 2 Data_2 : Minute, Data_3 : Hour (see Table B)

Table A: Day + Month + Year

Table B: Minute and Hour MSB LSB

Data_1 Data_0

Year (7-bit) Month (4-bit) Day (5-bit)

6 5 4 3 2 1 0 3 2 1 0 4 3 2 1 0

MSB LSB

Data_3 Data_2

Document Name: BG Meter multiple parameters Update: 2018/02/23

Document Version:1. 8

Proprietary and confidential. 8 / 16

2.2.2. [0x24] Read device model

Message Name: Read device model

Message ID: 0x24

Message Description: Read device model

Src / Des: GW / MD

Length: Request: 8 bytes / Response: 8 bytes

Format

Message name Direction CMD Data_0 Data_1 Data_2 Data_3

Read device model. GW MD 0x24 0x0 0x0 0x0 0x0

Response from device. GW MD 0x24 device model

Field Description

Field Size

(Bytes) Value Description.

CMD 1 0x24 GW commands MD to response device model

DATA 4 0x0

Field Size

(Bytes) Value Description.

ACK 1 0x24 MD response device model

DATA 2 Data_1 + Data_0 : Device model in four digits (Data_1 is the higher byte)

DATA 2 Data_3, Data_2 : undefined

2.2.3. [0x25] Read the storage data with index, part 1(time)

Message Name: Read the storage data with index, part 1(time) NA Hour (5 -bit) NA Minute (6 -bit)
4 3 2 1 0 5 4 3 2 1 0

Document Name: BG Meter multiple parameters Update: 2018/02/23

Document Version:1. 8

Proprietary and confidential. 9 / 16

Message ID: 0x25

Message Description: Read data (part 1: measure date and time)

Src / Des: GW / MD

Length: Request: 8 bytes / Response: 8 bytes

Format

Message name Direction CMD Data_0 Data_1 Data_2 Data_3

Read the storage data,

part1 . GW MD 0x25 Index 0x0 0x0

Response from device. GW MD 0x25 M_Date M_Time

Field Description

Field Size

(Bytes) Value Description.

CMD 1 0x25 GW commands MD to response measure data (part 1)

Index 2 0x00 is the latest reading in the device

DATA 2 0x00 NA

Field Size

(Bytes) Value Description.

ACK 1 0x25 MD response measure date, time and glucose type

M_Date 2 The measure date (see Table A in 2.3.1 [0x23])

M_Time 2 The measure time (see Table B in 2.3.1 [0x23])

2.2.4. [0x26] Read the storage data with index, part 2(result)

Message Name: Read the storage data with index, part 2(result)

Message ID: 0x26

Message Description: Read data (part 2: measure result)

Src / Des: GW / MD

Length: Request: 8 bytes / Response: 8 bytes

Format

Message name Direction CMD Data_0 Data_1 Data_2 Data_3

Read the storage data, GW MD 0x2 6 Index 0x0 0x0
Document Name: BG Meter multiple parameters Update: 2018/02/23
Document Version:1. 8

Proprietary and confidential. 10 / 16
part2.

Response from device. GW MD 0x26 Value Type I / Type II

Field Description

Field Size

(Bytes) Value Description.

CMD 1 0x26 GW commands MD to response measure data (part 1)

Index 2 0x00 is the latest reading in the device

DATA 2 0x00 NA

Field Size

(Bytes) Value Description.

ACK 1 0x26 MD response measure date and time

Value 2 The measure result (see Table C below)

The unit is mg/dL (mg/dL = mmol/L * 18)

Type I /

Type II 2 Parameters (see Table E below)

[Type I]

0x0: Gen

0x1: AC (before meal)

0x2: PC (after meal)

0x3: QC (quality control)

[Type II]

0x0: Gen.

Table C: Reading

Table D: Parameters MSB LSB

Data_1 Data_0

Value

15 14 13 12 11 10 9 8 7 6 5 4 3 2 1 0

MSB LSB

Data_3 Data_2

Type I Type II (4 bits) Code (10 bits)

1 0 3 2 1 0 9 8 7 6 5 4 3 2 1 0

Document Name: BG Meter multiple parameters Update: 2018/02/23

Document Version:1. 8

Proprietary and confidential. 11 / 16

2.2.5. [0x27] Read device serial number, part 1

Message Name: Read the device serial number, part 1

Message ID: 0x27

Message Description: Read the device serial number, part 1. This command cooperates with 0x28 to complete the serial number.

Src / Des: GW / MD

Length: Request: 8 bytes / Response: 8 bytes

The Serial Number for the device is composed of 16 characters (0~9, A~F) : SN_7 + + SN_0
(example: 427211713000001D, the former 4 bytes come from 0x28 and the later come from 0x27)

Format

Message name Direction CMD Data_0 Data_1 Data_2 Data_3
Read device serial
number, part 1. GW MD 0x27 0x0 0x0 0x0 0x0
Response from device. GW MD 0x27 SN_0 SN_1 SN_2 SN_3

Field Description

Field Size

(Bytes) Value Description.

CMD 1 0x27 GW commands MD to response the device serial
number , part 1
DATA 4 0x0

Field Size

(Bytes) Value Description.

ACK 1 0x27 Read the device serial number, part 1
DATA 4 Data_0 ~ Data_3 : SN_0 ~ SN_3

2.2.6. [0x28] Read device serial number, part 2

Message Name: Read the device serial number, part 1

Document Name: BG Meter multiple parameters Update: 2018/02/23

Document Version:1. 8

Proprietary and confidential. 12 / 16

Message ID: 0x2 8

Message Description: Read the device serial number, part 2. This command cooperates
with 0x2 7 to complete the serial number.

Src / Des: GW / MD

Length: Request: 8 bytes / Response: 8 bytes

The Serial Number for the device is composed of 16 characters (0~9, A~F) : SN_7 + + SN_0
(example: 427211713000001D, the former 4 bytes come from 0x28 and the later come from 0x27)

Format

Message name Direction CMD Data_0 Data_1 Data_2 Data_3

Read device serial

number, part 2. GW MD 0x28 0x0 0x0 0x0 0x0

Response from device . GW MD 0x28 SN_4 SN_5 SN_6 SN_7

Field Description

Field Size

(Bytes) Value Description.

CMD 1 0x28 GW commands MD to response the device serial
number, part 2
DATA 4 0x0

Field Size

(Bytes) Value Description.

ACK 1 0x28 Read the device serial number, part 2
DATA 4 Data_0 ~ Data_3 : SN_4 ~ SN_7

2.2.7. [0x2B] Read storage number of data

Message Name: Read storage number of data

Message ID: 0x2 B

Message Description: Read data counts

Src / Des: GW / MD

Length: Request: 8 bytes / Response: 8 bytes

Document Name: BG Meter multiple parameters Update: 2018/02/23
Document Version:1. 8

Proprietary and confidential. 13 / 16

Format

Message name Direction CMD Data_0 Data_1 Data_2 Data_3
Read Storage Number
of Stored Readings GW MD 0x2B 0x0 0x0 0x0 0x0
Response from device GW MD 0x2B Storage Number 0x0 0x0

Field Description

Field Size

(Bytes) Value Description.

CMD 1 0x2B GW commands MD to response the Number of
Stored Readings
DATA 4 0x0

Field Size

(Bytes) Value Description.

ACK 1 0x2B MD response the Number of Stored Readings
DATA 2 Data_1+Data_0 = Storage Number (Word)

2.2.8. [0x33] Write System Clock time (4 bytes)

Message Name: Write device clock time (4 bytes)

Message ID: 0x33

Message Description: Write the device date and time (to minute)

Src / Des: GW / MD

Length: Request: 8 bytes / Response: 8 bytes

Please refer to 0x23 command section for the format of date and time.

Format

Message name Direction CMD Data_0 Data_1 Data_2 Data_3
Write System Clock
time. GW MD 0x33 Day+Month+Year Minute Hour
Document Name: BG Meter multiple parameters Update: 2018/02/23
Document Version:1. 8

Proprietary and confidential. 14 / 16

Response from device. GW MD 0x33 Day+Month+Year Minute Hour

Field Description

Field Size

(Bytes) Value Description.

CMD 1 0x33 GW commands MD to write System Clock time (4
bytes)
DATA 2 Data_1 + Data_0 : Day+Month+Year (word)
DATA 2 Data_2 : Minute, Data_3 : Hour

Field Size

(Bytes) Value Description.

ACK 1 0x33 MD write System Clock time (4 bytes) O.K.
DATA 2 Data_1 + Data_0 : Day+Month+Year (word)
DATA 2 Data_2 : Minute, Data_3 : Hour

2.2.9. [0x50] Turn off the device

Message Name: Turn off the device

Message ID: 0x50

Message Description: Turn off the device from G W

Src / Des: GW / MD

Length: Request: 8 bytes / Response: 8 bytes

Format

Message name Direction CMD Data_0 Data_1 Data_2 Data_3

Turn off the device GW MD 0x50 0x0 0x0 0x0 0x0

GW MD 0x50 0x0 0x0 0x0 0x0

Field Description

Field Size

(Bytes) Value Description.

CMD 1 0x50 GW commands MD to turn off the power.

Document Name: BG Meter multiple parameters Update: 2018/02/23

Document Version:1. 8

Proprietary and confidential. 15 / 16

DATA 4 0x0 NA

Field Size

(Bytes) Value Description.

ACK 1 0x50 MD responses.

DATA 4 0x0 NA

2.2.10. [0x52] Clear/Delete all Memory

Message Name: Clear/Delete all memory

Message ID: 0x52

Message Description: Delete all measurement data in the device

Src / Des: GW / MD

Length: Request: 8 bytes / Response: 8 bytes

Format

Message name Direction CMD Data_0 Data_1 Data_2 Data_3

Clear/Delete all Memory

Storages GW MD 0x52 0x0 0x0 0x0 0x0

Response from device. GW MD 0x52 0x0 0x0 0x0 0x0

Field Description

Field Size

(Bytes) Value Description.

CMD 1 0x52 GW commands MD to Clear/Delete all Memory

Storages in Device EEPROM

DATA 4 0x0 NA

Field Size

(Bytes) Value Description.

ACK 1 0x52 Clear/Delete all Memory Storages in Device EEP ROM

O.K.

DATA 4 0x0 NA

Document Name: BG Meter multiple parameters Update: 2018/02/23

Document Version:1. 8

2.2.11. [0x54] Notification for entering communication mode

Message Name: Notification for entering communication mode

Message ID: 0x54

Message Description: Some models will send 0x54 command to host as a notification to indicate meter is ready for command from host. Host has to deal with this kind of notification in the beginning of connection (like RS-232) and prevent it to be mixed with the expected command response in receiving queue of host. This command only supports cable transmission, so Bluetooth transmission is not supported.

Src / Des : MD / GW

Length: Request: 8 bytes / Response: NA

Format

Message name	Direction	CMD	Data_0	Data_1	Data_2	Data_3
--------------	-----------	-----	--------	--------	--------	--------

Notification of entering						
--------------------------	--	--	--	--	--	--

communication mode	GW	MD	0x54	0x0	0x0	0x0
--------------------	----	----	------	-----	-----	-----

No need to response from GW side

Field Description

#	Field	Size
---	-------	------

(Bytes)	Value	Description.
---------	-------	--------------

CMD	1	0x54	MD reports to GW that it enters communication mode.
-----	---	------	---

DATA	4	0x0	NA
------	---	-----	----